



DataSheet

CATALOGUE #: 8FFS

PRODUCT NAME: Fatty acid binding protein (FABP) free serum

Source: Pooled normal human serum. Blood samples from donors were tested and found to be negative for HBsAg, HIV-1 and HIV-2 antibodies, and HCV.

Application: FABP free serum can be used for preparation of human FABP standards and calibrators.

Analysis: According to ELISA testing FABP free serum does not contain more than 0.5 ng/ml of human FABP.

Purification: Immunoaffinity chromatography.
Gel for affinity column was prepared by coupling anti-FABP monoclonal antibodies (several monoclonal antibodies specific for different epitopes of FABP molecule) with BrCN-activated Sepharose CL 4B.

Presentation: Frozen liquid.

Storage: -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi