



DataSheet

CATALOGUE #: 8P97

PRODUCT NAME: PAPP-A, homodimeric form (dPAPP-A), recombinant

Source: Expressed in human cell line

Application: Calibrator or standard in immunoassays, immunogen for antibody production.

Analysis: Purity > 90 % (SDS-PAGE in reducing conditions)
Immunoreactivity is checked using monoclonal antibodies specific to PAPP-A subunit.

Purification: Immunoaffinity purification.

Presentation: Frozen in 30 mM Tris-HCl, 90 mM NaCl, 40% Propandiol, pH 7.5

Storage: -70 °C (-65 ... -80 °C allowed)

Other information: Avoid repeated freezing and thawing.
It is recommended to aliquot the product.
This product is available only in research amounts.

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6
FI-20520 Turku FINLAND
www.hytest.fi | hytest@hytest.fi