



# DataSheet

**CATALOGUE #:** 4IG8

**PRODUCT NAME:** Monoclonal mouse anti-insulin-like growth factor binding protein-1 (IGFBP-1)

---

**MAbs:** G5F8, C7B9

**Specificity:** MAbs do not cross-react with human insulin or proinsulin

**MAb isotypes:** **IgG1** for MAb C7B9  
**IgG2a** for MAb G5F8

**Applications:** Antibodies are working in ELISA and Western blotting.  
Recommended pair for sandwich ELISA (capture-detection):  
C7B9 – G5F8

**Purification:** Salt fractionation and ion-exchange chromatography

**Presentation:** PBS, pH 7.4, 0.04 % sodium azide (NaN<sub>3</sub>)

**Storage:** +4 °C (+2 ... +8 °C allowed)

**Material safety note:** This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

**HyTest Ltd.**

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)