



# DataSheet

**CATALOGUE #:** 8IGB1

**PRODUCT NAME:** Insulin-like growth factor binding protein-1 (IGFBP-1), semipurified

---

**Source:** Pooled amniotic fluid  
Blood sample from the tissue donor was tested and found negative for HBsAg, HIV type 1 and 2 antibodies, HCV, and syphilis.

**Applications:** Calibrator or standard in immunoassays.

**Immunological activity:** Immunological activity of IGFBP-1 has been determined by IGFBP-1 IEMA test (Medix Biochemica, Finland).

**Purification:** Salt- fractionation and ion-exchange chromatography.

**Presentation:** Lyophilized, essentially salt free.

**Storage:** -20 °C (-15 ... -30 °C allowed)

**Material safety note:** This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

**HyTest Ltd.**

Intelligate 1, 6th floor, Joukahaisenkatu 6  
FI-20520 Turku FINLAND  
[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)