



# DataSheet

**CATALOGUE #:** 8CA19

**PRODUCT NAME:** CA19-9

---

**Source:** Human metastatic liver carcinoma.  
Blood sample from the tissue donor was tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies, HCV, and syphilis.

**Applications:** Immunogen, calibrator in immunoassays, in vitro studies.

**Specific activity:** Specific activity confirmed by enzyme immunoassay.

**Presentation:** PBS, pH 7.4, 3 % sucrose, 0.05 % of sodium azide (NaN<sub>3</sub>).

**Storage:** -20 °C (-15 ... -30 °C allowed)

**Material safety note:** This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.  
Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

**HyTest Ltd.**

Intelligate 1, 6th floor, Joukahaisenkatu 6  
FI-20520 Turku FINLAND  
[www.hytest.fi](http://www.hytest.fi) | [hytest@hytest.fi](mailto:hytest@hytest.fi)