



DataSheet

CATALOGUE #: 8C29

PRODUCT NAME: CA-125

Source: Human adenocarcinoma.
Blood sample from the tissue donor was tested and found negative for HBsAg, HIV-1 and HIV-2 antibodies, HCV, and syphilis.

Application: Immunogen, in vitro studies.

Immunological specificity: Immunological specificity confirmed by reaction with monoclonal antibody that is specific for CA-125.

Presentation: 20 mM Tris-HCl, pH 7.5, 0.1 % NaCl, 0.09 % NaN₃

Storage: -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.
Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6
FI-20520 Turku FINLAND
www.hytest.fi | hytest@hytest.fi