



DataSheet

CATALOGUE #: 4SA11

PRODUCT NAME: Monoclonal mouse anti-serum amyloid A (SAA)

MAbs: SAA1, SAA6, SAA15

Hybridoma clones have been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with human recombinant SAA, conjugated with carrier protein.

Specificity: Human recombinant SAA

MAb isotypes: IgG1 for MAbs SAA1, SAA6 and SAA15

Applications: SAA immunoassay. All MAbs recognize SAA in ELISA and Western blotting.

Recommended pairs for sandwich immunoassay (capture-detection):

SAA1 - SAA15

SAA6 - SAA15

SAA15 - SAA1

SAA6 - SAA1

SAA1 - SAA6

Purification: Chromatography on protein A Sepharose

Presentation: PBS, pH 7.4, 0.09 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi