



DataSheet

CATALOGUE #: 8CFS

PRODUCT NAME: CRP free serum

Source: Pooled normal human serum. Blood samples from donors were tested and found to be negative for HBsAg, HIV-1 and HIV-2 antibodies, and HCV.

Application: CRP free serum can be used for preparation of CRP standards and calibrators.

Analysis: According to ELISA testing CRP free serum does not contain more than 0.02 µg/ml of human CRP.

Purification: Immunoaffinity chromatography of pooled normal human serum.
Gel for affinity column was prepared by coupling anti-CRP monoclonal antibodies with BrCN-activated Sepharose CL 4B.

Presentation: Frozen liquid.

Storage: -20 °C (-15 ... -30 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.