



DataSheet

CATALOGUE #: RIF17

PRODUCT NAME: Monoclonal mouse anti-influenza virus B group antigen

MAbs: R2/3

MAb isotypes: IgG2a for MAb R2/3

Applications: Indirect EIA, Western blotting

Purification: Chromatography on protein A Sepharose

Presentation: PBS, pH 7.4, 0.09 % sodium azide (NaN₃)

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.

HyTest Ltd.

Intelligate 1, 6th floor, Joukahaisenkatu 6

FI-20520 Turku FINLAND

www.hytest.fi | hytest@hytest.fi