



DataSheet

CATALOGUE #: 8HS7-3

PRODUCT NAME: HBsAg, recombinant, adr subtype

Source: Yeast

Presentation: 10 mM phosphate, 145 mM NaCl , 0.09 % sodium azide (NaN₃). The product is sterile filtered.

Analysis: Purity > 98 % (HPLC)

Immunological activity: Tested in ELISA with anti-HBsAg antibodies.

Application: ELISA (coating and conjugate)
Immunochromatography (capture and conjugate)
Immunogen for monoclonal antibody production

Storage: +4 °C (+2 ... +8 °C allowed)

Material safety note: This product is sold **for research use only**. Standard Laboratory Practices should be followed when handling this material.

Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.