




DataSheet

CATALOGUE #: 4CA19

PRODUCT NAME: Monoclonal mouse anti-CA19-9

MAbs:	X91 Hybridoma clone has been derived from hybridization of Sp2/0 myeloma cells with spleen cells of Balb/c mice immunized with affinity purified human CA19-9 antigen.
Specificity:	CA19-9 antigen
MAb isotypes:	IgM for MAb X91
Applications:	Detection of CA19-9 antigen in immunoassays. MAbs are working in immunohistochemistry (see picture below).
Purification:	Ion exchange chromatography
Presentation:	0.15 M NaCl, 10 mM Tris, pH 7.5, 0.05 % sodium azide (NaN ₃)
Storage:	+4 °C (+2 ... +8 °C allowed)
Material safety note:	This product is sold for research use only . Standard Laboratory Practices should be followed when handling this material. Product contains sodium azide as a preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.



MAb X91: Immunostaining of CA19-9 antigen in normal pancreas